

TRANSFER AND EMERGENCY STRETCHERS WHEEL CHAIR



AGENCINOX
LA SANTÉ EXIGE LA QUALITÉ

AXIOME

aviso

ABSOLUTE

AMBULYS

CE

Variable height



AXIOME

Variable height	Trendelenburg
Width 500mm	Width 500mm
CB501	CB502
CB501MC	CB502MC
CB501 5R	CB502 5R
CB501MC5R	CB502MC5R
Width 600mm	Largeur 600mm
CB601	CB602
CB601MC	CB602MC
CB601 5R	CB602 5R
CB601MC5R	CB602MC5R
Width 700mm	Width 700mm
CB701	CB702
CB701MC	CB702MC
CB701 5R	CB702 5R
CB701MC5R	CB702MC5R
Width 800mm	Width 800mm
CB801	CB802
CB801MC	CB802MC
CB801 5R	CB802 5R
CB801MC5R	CB802MC5R

- Height adjustable by hydraulic spring, left and right operated
- Removable mattress. Flame retardant M1 (83% PVC / 17% PES). Resin hard top
- Collapsible side rail
- Adjustable backrest (from 0 to 90°). Double gaz spring assisted.
- Telescopic and collapsible serum rod – 2 hooks.
- Accessories holders and bumpers.
- 4 wheels 200 mm diam. Optional 5 th wheel for better manoeuvrability
- 2 retractables handles
- Epoxy (grey 7035 and 9007). Corrosion resistant – tropicalized
- Optional trendelenburg position by double gaz spring . “Easy lift system”

Trendelenburg



Knee section




Knee section
Mattress thickness 120mm
Width 600mm
CB602P
CB602P 5R
Width 700mm
CB702P
CB702P 5R

Accessories

AXIOME

- Ref. CBSO Oxygen holder
- Ref. CBPA Basket
- Ref. CBSS Monitoring holder with straps
- Ref. CBSR Paper roll holder
- Ref. CBSU Urinal bag holder
- Ref. CBTS Additionnal telescopic serum rod
- Ref. CBTSR2 Additionnal telescopic foldable serum rod (head end only)
- Ref. CBTR Food tray
- Ref. CBPB Side rails protection
- Ref. CBPR File holder
- Ref. CBK7 Radio holder including the full length Translatory motion
- Ref. CBSM Over bed radiotransluscent mattress
- Ref. CBARE Shoulders restraining strap
- Ref. CBARM Ankle restraining strap
- Ref. CBARP Wrist restraining strap
- Ref. CBPB Side rail protection
- Ref. CBT Ophtalmic surgery headrest
- Ref. CBPP Additional pushing handle



Technical specifications

Mattress width	Mattress length	Mattress thickness	Mattress Thickness (knee sect.)	Overall width	Overall length	Low position	Hight position	Weight without accessories	Load capacity
480mm	1860mm	80mm		614mm	2040mm	580mm	890mm	70kg	250kg
580mm	1860mm	80mm	120mm	714mm	2040mm	580mm	890mm	80kg	250kg
680mm	1860mm	100mm	120mm	814mm	2040mm	580mm	890mm	100kg	250kg
780mm	1860mm	100mm		914mm	2040mm	580mm	890mm	110kg	300kg



Model presented : Aviso Ref. CBL601MC with CBLPPE



Model presented : Aviso Ref. CBL6025R with CBLTS and CBLPPE with trendelenburg

Configuration and specifications

- Height adjustable by hydraulic spring.
- Trendelenburg version available.
- Removable mattress, flame retardant M1, 35 kg/m³ housse 83% PVC / 17% PES, resin hard top.
- Collapsible side rail.
- Adjustable backrest (from 0 to 90°). Double gaz spring assisted (easy lift).
- Accessories holders and bumpers.
- 4 wheels Ø150 mm. Central braking system left/right side operated. (optional wheels Ø200mm)
- Option : 5th wheel for better manoeuvrability
- Epoxy coated frame (grey 7035 and 9007). Corrosion resistant - tropicalized.
- 200 kg capacity.

Transfer mattress available, PVC coated coton cover.

Optional accessories

- Basket / Shelf **Ref. CBLPA**
- 2 retractables handles **Ref. CBLPPE**
(head or foot end)
- 2 fixed handles **Ref. CBLPPF**
(head or foot end)
- Removable serum rod **Ref. CBLTS**
- Telescopic and foldable serum rod **Ref. CBLTSR2**
(head end only)
- Removable telescopic serum rod **Ref. CBLTSA**
- Holder for one O2 **Ref. CBSO**
- Urinal bag holder **Ref. CBLSU**
- Chart/radio holder **Ref. CBLPR**
- Paper roll holder **Ref. CBLSR**
- 2 aluminium rail **Ref. CBLRA**
- Clamp **Ref. CBLEF**

Dimensions

Height	Trendelenburg
Width 600mm	Width 600mm
CBL601	CBL602
CBL601MC	CBL602MC
CBL601 5R	CBL602 5R
CBL601MC5R	CBL602MC5R
Width 700mm	Width 700mm
CBL701	CBL702
CBL701MC	CBL702MC
CBL701 5R	CBL702 5R
CBL701MC5R	CBL702MC5R

Technical specifications

Mattress width	Mattress length	Mattress thickness	Overall width	Overall length	Low position	High position	Weight without accessories
580mm	1860mm	80mm	717mm	2040mm	525mm	917mm	80kg
680mm	1860mm	100mm	817mm	2040mm	525mm	917mm	95kg



Comfortable pushing handle

Electrical backrest

Control pad 4 adjustments



Shelf utility tray

Ø 100 wheel with brake



Electrical leg-rest
Or shortened leg-rest



FA.GF
Grey / framboise
FAR.GF



FA.LP
Lavande / prune
FAR.LP



FA.GO
Grey / orange
FAR.GO



FA.AM
Anis / brown
FAR.AM

Configurations and specifications

Patient chair, electrical with 4 « LINAK » actuators allowing height adjustment from 560 mm to 920 mm, Trendelenburg / reverse 20°/18°, backrest and legrest from 90° to 180°. All adjustments are done by hand control pad.

Optionnal shortened leg-rest (373mm).

Features:

- Basket for accessories.
- 4 x 100mm wheels including 2 with brake. (non marking).
- Epoxy grey coated 7035 RAL , frame finishing, anti rust treatment, corrosion resist tropicalized.
- Max. load 200 kg.
- Saddle : choice of 4 bi-colors upholstery.
- Plastic cover for leg-rest

3 parts couch with additional memory foam, 100mm thick 35 kg/m³, PVC coating upholstery: Abrasive resistant testing Martindale 120 000 cycles. M2 Fire proof form.

Clean regularly with a damp soapy cloth. Do not use strong detergents, solvents cleaners, polish, wax or patent chemicals: harmful for the material. Some stains like ball pen, food coloring should be cleaned off immediately to avoid any permanent absorption into the upholstery.

Optionnal accessories (*assembly plant*)

- | | |
|---|---------------------|
| • Revolving arm-rest | Ref. FAAR |
| • Adjustable blood taking arm-rest on ball | Ref. FAABR |
| • 2 stainless steel rail 25x10mm | Ref. FARA |
| • Paper roll holder | Ref. FASR |
| • Urinal bag holder | Ref. FASU |
| • Battery | Ref. FABAT |
| • Stainless steel foot-rest | Ref. FARP |
| • 4X125mm directional and central braking
Wheels | Ref. FAC+color ref |
| • Shortened leg-rest | Ref. FARC+color ref |
| • Additional Brake on front wheels | Ref. FAPF |

Dimensions

Couch :

- Length 1900mm
- Width 550mm

Ext. dimensions:

- Length 1900mm
- Width 820mm
(*arm-rest version*)

Height:

- Position low 560mm
- Position high 920mm

Gross weight : 80 kg

Other optionnal accessories

- | | |
|-------------------------------|-----------|
| • Armrest clamp Ø 16mm | Ref. FAEF |
| • Multi use shelf | Ref. FAPC |
| • Serum rod 2 hooks | Ref. FATS |
| • Flexible control pad holder | Ref. FAST |
| • Adjustable head pillow | Ref. FART |
| • LED reading light | Ref. FALL |

AMBULYS

MULTI PURPOSE CHAIR
Mattress size 600 or 700mm



Pushing handles,
 Seat or flat
 position



Safety side rail
 extendable 1000 to
 1150mm (option)

Control pad 4
 adjustments

Electrical
 back-rest

Electrical
 leg-rest

Emergency
 position control
 device

Ø125 wheel or Ø150
 (with brake)



FHJ.GF
 Grey /framboise



FHJ.LP
 Lavande /prune



FHJ.GO
 Grey /orange



FHJ.AM
 Anis /brown

Configuration and specifications

- 4 electrical adjustments with control pad : Height, back-rest, leg-rest, trendelenburg, reverse trendelenburg
- 2 Pushing handles for seat and/or flat position.
- Removable safety side rails (optional), length 1000mm or extendable side rails from 1000 to 1150mm, in accordance with the standard bed-ins, useful for adults and children.
- Emergency position control.
- Accessories holder.
- 4 wheels with brake Ø 125mm 2 side control (non marking).
- Steel frame, Epoxy grey coated 7035, corrosion resist tropicalized.
- 200 kg loading capacity.
- Saddle: choice of 4 bi-colors upholstery.

Saddlery with memory foam depth 100mm 35 kg/m³, PVC coating upholstery: Abrasive resistant testing Martindale 120 000 cycles. M2 Fire proof.

Clean regularly with a damp soapy cloth. Do not use strong detergents, solvents cleaners, polish, wax or patent chemicals: harmful for the material. Some stains like ball pen, food colorings should be cleaned off immediately to avoid any permanent absorption into the vinyl.

Optionnal accessories (*assembly plant*)

- | | |
|---|--------------|
| • 2 rails inox 25x10mm | Ref. FHJRA |
| • 4 wheels Ø150mm central braking | Ref. FHJR150 |
| • Urinal bag holder | Ref. FHJSU |
| • Battery | Ref. FHJBAT |
| • Stainless steel foot-rest | Ref. FHJCPF |
| • Stainless steel adjustable foot-rest | Ref. FHJCPR |
| • Safety side rails Length 1000 mm | Ref. FHJBF |
| • Safety side rails, extendable 1000-1150 | Ref. FHJBE |

Specifications / Dimensions


	FHJ600	FHJ700
Mattress		
Length	2100	2100
Width	600	700
Overall		
Length	2190	2190
Width	780	880
height		
Low position	540	540
High position	905	905
Weight	90kg	96kg

Other optionnal accessories

- | | |
|-------------------------------|------------|
| • Clamp Ø 16mm | Ref. FHJEF |
| • Multi use shelf | Ref. FHJPC |
| • Serum rod 2 hooks | Ref. FHJTS |
| • Flexible control pad holder | Ref. FHJST |
| • Adjustable head pillow | Ref. FHJRT |
| • LED reading light | Ref. FHJLL |
| • 2 armrest mounted on ball | Ref. FHJAB |
| • Brake on front wheel side | Ref. FHJPF |

ABSOLUTE **AMBULYS**


Comfortable - Safety - Ergonomic - Handling hability - Reliability - Sturdiness



WHEEL CHAIR OPTIONS



FIND OTHER PRODUCTS AGENCINOX:





For more information or to order, call



our Sole Distributor in Ireland

Pat Coakley Medical Ltd

Units E3 E4
Network Enterprise Park
Kilcoole A63 H244
Co Wicklow
IRELAND

+353 1 2018070

sales@patcoakleymedical.ie



AGENCINOX

Z.I. 89110 Aillant sur Tholon – France
Export : Tél : 33 3 86 92 39 27. Fax 33 3 86 91 54 42

<http://www.agencinox.com>

e mail : contact@agencinox.com

S.A.S au capital de 800 000 € - RC Sens B 353 427 495
Siret 353 427 495 00018 – NAF 3250A

