

Stainless Steel Dins



Stainless Steel Dins

Our stainless steel range of DINS are a quality product available as a Flat Base Basket with Perforated Sides or Full Wire Mesh Baskets. The Flat Base with Perforated Sides Basket prevents finer, delicate instruments from slipping through. Both versions have 2 inner mounted drop down handles.

- > Manufacture to DIN 58952
- > Electro-Polished
- > Spot-welded at cross points
- > 304 Stainless Steel
- > Burr Free
- > Cross-hatched flat mesh base



Flat Base Baskets with Perforated Sides suffix *PER

Full Wire Mesh Baskets suffix *MESH

Please add the suffix *PER or *MESH when you place your order to identify the tray you require

Ref	Dimensions (LxWxD mm)	Ref	Dimensions (LxWxD mm)
1001	240 x 250 x 60	1004	410 x 250 x 60
1002	240 x 250 x 75	1005	410 x 250 x 75
1003	240 x 250 x 100	1006	410 x 250 x 100

Ref	Dimensions (LxWxD mm)	Ref	Dimensions (LxWxD mm)
1007	480 x 250 x 60	1010	540 x 250 x 60
1008	480 x 250 x 75	1011	540 x 250 x 75
1009	480 x 250 x 100	1012	540 x 250 x 100



Pat Coakley Medical Ltd.

Units E3 E4 Network Enterprise Park
Kilcoole A63 H244 Co Wicklow
IRELAND

+353 12018070

sales@patcoakleymedical.ie

**Pat Coakley
Medical**

