



**NEOTECH®**

## ***Making a Difference***



# Product Catalog

neotechproducts.com

# New Product Development

Do you have an idea for a new product? Neotech can help make it a reality!  
Do you have a problem in your unit, but no idea what the product solution is?  
We can help with that too!

We have a great history of successfully licensing new product ideas and bringing them to market. We can help you every step of the way. We're easy to work with and we move forward quickly at no cost to you. Visit [neotechproducts.com/newproduct](http://neotechproducts.com/newproduct) to learn more or call us to receive a new product submission form.



DO YOU HAVE A PRODUCT IDEA?

NEOTECH CAN HELP!



Look for this icon for our Home Health products.  
Visit [neotechhomehealth.com](http://neotechhomehealth.com) for more information.

NeoHug®	2	Suction Caddy™	19
NeoHeart™	3	NeoShades® with Tabs	20
NeoLead®	4	NeoShades® with Headstrap	21
Radiolucent NeoLead®	5	EZ-Hold®	22
NeoGlo® & NeoGlo Sleeve	6-7	NeoBond®	23
Neotech RAM Cannula®	8-9	Hold-A-Line™	24
EZCare Softouch™	10	NeoBridge®	25
Neotech ChinStrap™	11	NeoSmile™	26
NeoPulse™	12	SoftSil™	27
NeoBar®	13	NeoPockets™	28
NeoSeal®	14	NeoSpoon™	29
Neotech Meconium Aspirator®	15	Neotech Cell Wipe™	30
Little Sucker®	16-17	Neotech Cell Shield™	31
Little Sucker Cover®	16-17	NeoGrip®	32
Little Sucker® Nasal Tip	16-17	Clinical Training and Research	33
NeoSucker® Curved	18		

# About NeoBond® Hydrocolloid

## Important, Please Read!

Skin friendly NeoBond Hydrocolloid is a protective adhesive that helps decrease the risk of damage to the fragile epidermis of the neonate. AWHONN Neonatal Skin Care (Third Edition, 2013) recommends that hydrocolloids be used in place of other adhesives whenever possible.

Heat is the catalyst for hydrocolloid. The following steps are necessary to ensure the effectiveness of hydrocolloid based products.



### PREPARATION

Clean and dry skin per hospital protocol.

Note: It is important for skin to be as clean and dry as possible. Products will not adhere properly to moist skin or hair. Oils and lotions will also adversely affect adhesion.

**DO NOT USE ALCOHOL.**



### APPLICATION

Before removing liners, warm hydrocolloid in hands for 60 seconds. May also be held under radiant or heel warmer for 10 to 15 seconds.

Remove liners, apply hydrocolloid base to skin and hold in place for 60 seconds to ensure proper adhesion.



### REMOVAL

Replace product every 5 days or per hospital protocol, whichever is sooner.

Saturate hydrocolloid with water or saline.

Slowly peel away from skin as you swab with water or saline.



# NeoHug®

## Utility Device Holder



N680BL-S



N680BO-S



N680BP-S



N680GR-S



N680VL-S

### FEATURES AND BENEFITS

- Easy to use and attaches to most smooth surfaces
- Large suction cup for strong hold and easy removal
- Holds under heat and humidity
- Reposition repeatedly with no adhesive residue
- Bendable arms to hold a variety of light weight items
- Ideal for holding pacifier, suction tip, tape and more
- Playful chimp design available in five bright colors
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Color	Qty/Unit
N680BL-S	Blue	25/box
N680BO-S	Bright Orange	25/box
N680BP-S	Bright Pink	25/box
N680GR-S	Green	25/box
N680VL-S	Violet	25/box
N680A-S	Assorted (5 of each color)	25/box


Patent Pending

# NeoHeart™

## ECG Pad



N335




N336

### FEATURES AND BENEFITS

- 3-in-1 pre-wired electrode array
- Pre-measured sensor position
- Radiolucent material nearly invisible on x-rays
- One touch application ideal for use in delivery room
- Larger surface area for enhanced adherence
- Engineered for back application, allowing chest to remain exposed
- Multi-colored wires for easier identification
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Size	Qty/Unit
N335	Small	20 sets/box
N336	Large	20 sets/box



Made in USA Patents Pending




# NeoLead®

## ECG Electrodes



N301



N300

### FEATURES AND BENEFITS

- Pre-wired electrodes
- Hydrogel center
- NeoBond® Hydrocolloid border designed to:
  - prevent hydrogel from oozing or drying prematurely
  - provide sensor stability for better tracings
  - minimize slipping
- Multi-colored wires for easier identification
- Not made with natural rubber latex or plasticizer DEHP
- Sets of three individually packaged

Cat No.	Item	Size	Qty/Unit
N300	Mini	3/4" circle	20 sets/box
N301	Standard	1" Circle	20 sets/box



Made in USA

U.S. Patent #6,795,722

# Radiolucent NeoLead®



Clinician Invented

## ECG Electrodes



N305



N306

### FEATURES AND BENEFITS

- Pre-wired electrodes are nearly invisible on x-rays
- Hydrogel center
- NeoBond® Hydrocolloid border designed to:
  - prevent hydrogel from oozing or drying prematurely
  - provide sensor stability for better tracings
  - minimize slipping
- Multi-colored wires for easier identification
- N305 set of three leads
- N315 set of five leads for two channel, continuous monitoring
- Not made with natural rubber latex or plasticizer DEHP
- Sets of leads individually packaged

### Micro NeoLead®

- Pre-wired electrodes are nearly invisible on x-rays
- Hydrogel center
- NeoBond® Hydrocolloid ends designed to:
  - provide sensor stability for better tracings
  - minimize slipping
  - reduce need for band or strap
- Multi-colored wires for easier identification
- Not made with natural rubber latex or plasticizer DEHP
- Sets of three individually packaged

Cat No.	Item	Size	Qty/Unit
N305	Radiolucent	1" Circle	20 sets/box
N306	Micro	1" x 3/8"	20 sets/box
N315	Pediatric	1" Circle	20 sets/box



Made in USA

U.S. Patent #6,795,722; #8,798,709

# NeoGlo®

Transilluminator







## FEATURES AND BENEFITS

- LED lights cool to the touch for patient safety
- Multiple light settings for user preference
- Ergonomic design for clinician comfort
- Compact size for use in incubator
- Illuminates veins, arteries and other internal structures
- Engineered to shut off before lights become ineffective
- Requires single AA alkaline battery (included)
- Eight colors to choose from
- Lanyard hoop
- Affordably priced
- Not made with natural rubber latex or plasticizer DEHP

### NeoGlo Sleeve Also Available



### NeoGlo Sleeve

- Protective cover helps to keep NeoGlo clean and protect against bodily fluids
- Durable and disposable
- Self-sealing flap
- Not made with natural rubber latex or plasticizer DEHP

Cat No.	Item	Qty/Unit
N555BL	Blue	1/box
N555AQ	Aqua	1/box
N555VL	Violet	1/box
N555BP	Bright Pink	1/box
N555PK	Rose	1/box
N555WH	White	1/box
N555SL	Silver	1/box
N555BK	Black	1/box
N55100	Sleeve 100Pack	100/box*


\*NeoGlo Sleeve: minimum purchase of 10 boxes



Made in USA Patent Pending

# Neotech RAM Cannula®

Nasal Oxygen Cannula





## FEATURES AND BENEFITS

- For use with low/high flow humidified oxygen
- Soft, curved prongs
- Color coded sizes
- 15 mm oxygen tubing adapter included
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Color	HCPCS	Qty/Unit
N4900	White	A4615	10/box
N4901	Green	A4615	10/box
N4902	Blue	A4615	10/box
N4903	Orange	A4615	10/box
N4904	Yellow	A4615	10/box
N4905	Teal	A4615	10/box
N4906	Purple	A4615	10/box
N4820	Oxygen Tubing Adapter		50/box



U.S. Patent #8,353,296; #D681,193;  
#9,027,562; #9,808,593; #9,861,777  
Additional Foreign Patents Pending

Featuring NeoFoam® Ultra

# EZCare™ Softouch


Tracheostomy Tube Holder



Parent Invented



Home Health



N926 18" 2-piece can be cut to size

## FEATURES AND BENEFITS

- Thick NeoFoam® Ultra for better stability
- Soft, breathable, skin friendly NeoFoam wicks perspiration
- Designed for comfort
- Lightweight and form fitting
- N926 18" can be cut to size
- Unique material helps to prevent sliding
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Size	HCPCS	Qty/Unit
N9105	5" Single piece	A7526	20/box
N9106	6" Single piece	A7526	20/box
N9107	7" Single piece	A7526	20/box
N9108	8" Single piece	A7526	20/box
N9109	9" Single piece	A7526	20/box
N9110	10" Single piece	A7526	20/box
N9111	11" Single piece	A7526	20/box
N925*	11" 2-piece (tan only)	A7526	20/box
N926**	18" 2-piece (tan only)	A7526	20/box
N926BK-10**	18" 2-piece (black only)		10/pack

Tan=TN, Camouflage=CO, Pink=PK, Blue=BL

\*1/2" tall

\*\*1" tall



Made in USA

# Neotech ChinStrap™

## Respiratory Support Accessory



N784

### FEATURES AND BENEFITS

- Thick, soft, breathable, skin friendly NeoFoam® wicks perspiration
- Flexibility enables minimal mouth movement
- Form fitting and adjustable
- Split chin design for comfort and stability
- Ear openings to prevent folding over
- Designed to reduce the risk of pressure ulcers or necrosis
- Unique material helps to prevent sliding
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged



Cat No.	Size	HCPCS	Qty/Unit
N784	Nano	A7036	20/box
N785	Micro	A7036	20/box
N786	Small	A7036	20/box
N787	Large	A7036	20/box
N788	XL	A7036	20/box
N789	Macro	A7036	20/box
N790	Jumbo	A7036	20/box



Made in USA U.S. Patent #10,342,693

# NeoPulse™

## Disposable Pulse Oximeter Wrap



Clinician Invented



HomeHealth



N781



N782



N783



### FEATURES AND BENEFITS

- Soft, breathable, skin friendly NeoFoam® wicks perspiration
- Lightweight, form fitting and adjustable
- Designed for comfort and stability
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Size	Dimensions	Qty/Unit
N781	Small	3/4" x 3"	100/box
N782	Large	1" x 5"	100/box
N783	XL	2" x 6"	100/box




Made in USA



# NeoBar®

## ET Tube Holder



### FEATURES AND BENEFITS

- Skin friendly NeoBond® Hydrocolloid tabs
- Macro and Jumbo sizes available with medical grade acrylic adhesive
- Designed to reduce extubations, help prevent palate trauma and allow for better oral care
- Eliminates need for tape near nose or mouth
- Can be cut for emergency removal
- Color coded sizes
- Not made with natural rubber latex or plasticizerDEHP
- Individually packaged

\*Improved curve design for N709 Mini, N710 Ultra and N711 Micro

Cat No.	Size	Color	Qty/Unit
N709	Mini	Purple	5/box
N710	Ultra	Yellow	5/box
N711	Micro	White	5/box
N712	Small	Green	5/box
N713	Large	Blue	5/box
N714	XL	Rose	5/box
N715H	Macro	Clear	5/box
N716H	Jumbo	Gray	5/box
With Medical Grade Acrylic Adhesive			
N715F	Macro	Peach	5/box
N716F	Jumbo	Carbon	5/box
Assorted			
N718 Transport (1 each N709,N710,N711, N712, N713, N714)			
N719N Assorted (2 each N709,N710,N711, N712, N713, N714)			





Made in USA      U.S. Patent #D434,496;#6,050,263

# NeoSeal<sup>®</sup>



## Nasal Prongs Seal




N420  
Actual Size



N421  
Actual Size



N422  
Actual Size



### FEATURES AND BENEFITS

- Designed to improve the seal for more effective respiratory support
- Protects the skin of the septum and nares
- Compatible with most nasal prongs
- Not made with natural rubber latex or plasticizer DEHP


Cat No.	Size	Qty/Unit
N420	Micro	150/box
N421	Small	150/box
N422	Large	150/box



Made in USA Patent Pending

# Neotech Meconium Aspirator®

## Meconium Removal Device




### FEATURES AND BENEFITS

- Designed to eliminate mouth suctioning of meconium
- Suction control port for intermittent suctioning
- Attaches to most endotracheal tubes
- Featured in the Textbook of Neonatal Resuscitation<sup>1</sup>
- N0102 (yellow) allows for easier identification
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

1) Textbook of Neonatal Resuscitation, 6th Edition (2011) pg. 43, Figure 2.7

Cat No.	Item	Qty/Unit
N0101	Clear	40/box
N0102	Yellow	40/box



*The original USA made*

# Little Sucker®

Oral and Nasal Suction Device



N204

# Little Sucker Cover®

Suction Device Protective Cover



N205C

# Little Sucker® Nasal Tip

Nasal Suction Device



N225



N205



N224

## FEATURES AND BENEFITS

### Little Sucker

- Designed for single handed suctioning
- Thumb port for intermittent suctioning
- Soft, flexible tip similar to a bulb syringe
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

### Little Sucker Cover

- Designed to protect the Little Sucker from damage, debris and external contaminants
- Bright blue color for better visibility
- Noise reducing port plug
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

### Little Sucker Nasal Tip

- Specialized tip designed to minimize potential tissue trauma
- Designed for single handed suctioning
- Thumb port for intermittent suctioning
- Soft, flexible tip similar to a bulb syringe
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Item	HCPCS	Qty/Unit
N204	Preemie	A4628	50/box
N204TP	Preemie 2-piece*	A4628	50/box
N204C	Preemie with Cover	A4628	50/box
N205	Standard	A4628	50/box
N205C	Standard with Cover	A4628	50/box
N205TP	Standard 2-piece*	A4628	50/box
N215	Cover	A4628	50/box
N224	Preemie Nasal Tip	A4628	50/box
N224TP	Preemie Nasal 2-Piece	A4628	50/box
N225	Standard Nasal Tip	A4628	50/box
N225TP	Standard Nasal 2-piece	A4628	50/box

\*2-piece Little Suckers are not compatible with the Cover



U.S. Patent #7,635,361;  
#D590,056; #D655,812; #D668,757;  
Additional Patents Pending

*Anatomically curved design*

# NeoSucker® Curved

Oral and Nasal Suction Device



## FEATURES AND BENEFITS

- Anatomically curved shape
- Designed for single handed suctioning
- Thumb port for intermittent suctioning
- Soft, flexible tip similar to a bulb syringe
- N208 XL ideal for older or larger pediatric patients
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Size	HCPCS	Qty/Unit
N206	Preemie	A4628	50/box
N207	Standard	A4628	50/box
N208	XL	A4628	50/box



U.S. Patent  
#7,635,361; #D590,056;  
#6,958,050; #D474,269



# Suction Caddy™


Suction Device Holder



## FEATURES AND BENEFITS

- Easy to use and attaches to most smooth surfaces
- Keeps suction devices and tubing close at hand and off the floor
- Convenient, designated write-on area for patient name and date
- Clear plastic allows for better visibility
- Tubing anchor helps to reduce noise
- Large suction cups provide strong hold and allow for easy removal
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Item	Qty/Unit
N245	Suction Caddy	20/box



Made in USA U.S. Patent #7,601,001

# NeoShades® with Tabs

Phototherapy Eye Shields



## FEATURES AND BENEFITS

- Soft, padded eye shield material with cute sunglasses design
- Ocular pockets designed for patient comfort and safety (N724, N725, N726)
- Skin friendly NeoBond® Hydrocolloid tabs hold NeoShades in place
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged



Cat No.	Item	Qty/Unit
N720	Micro	50/box
N721	Small	50/box
N722	Large	50/box
With Ocular Pockets		
N724	Micro	50/box
N725	Small	50/box
N726	Large	50/box



U.S. Patent #8,070,288;  
Additional Patents Pending

# NeoShades® with Headstrap

Phototherapy Eye Shields



## FEATURES AND BENEFITS

- Soft, padded eye shield material with cute sunglasses design
- Ocular pockets designed for patient comfort and safety
- Soft, patented NeoFoam® headstrap holds NeoShades in place
- Adjustable headstrap engineered to minimize slipping and help reduce head molding
- Unique headstrap shape designed to ensure ideal fit
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Item	Qty/Unit
N727	Micro	50/box
N728	Small	50/box
N729	Large	50/box



U.S. Patent #8,070,288;  
Additional Patents Pending




# EZ-Hold®


## Tubing and Cannula Holder



N756



N757



N758



### FEATURES AND BENEFITS

- Skin friendly NeoBond® Hydrocolloid base
- Fold over flap engineered to relieve stress on tubes
- Flap designed to be opened and closed repeatedly
- Clear flap allows for better visibility
- Mini size and shape ideal for micro preemies
- Not made with natural rubber latex or plasticizer DEHP

Cat No.	Item	HCPCS	Qty/Unit
N756	Mini	E1399	100/box
N757	Standard	E1399	100/box
N758	XL	E1399	100/box



Made in USA



U.S. Patent #D506,547

# NeoBond®

Hydrocolloid Adhesive Strips



N741



N742

## FEATURES AND BENEFITS


- Skin friendly alternative to acrylic tape
- Developed for long term attachment
- Ideal for use with limb boards
- Not made with natural rubber latex or plasticizer DEHP

Cat No.	Size	HCPCS	Qty/Unit
N741	1/2" x 6" Strips	A6234	100/box*
*50 sets of 2			
N742	1" x 15' Roll	A6234	1/box



# Hold-A-Line™

## Catheter and Tubing Holder



### FEATURES AND BENEFITS

- Skin friendly NeoBond® Hydrocolloid base
- Fold over flap designed to relieve stress on tubes and help prevent dislodgement
- Unique shape allows for protruding umbilical stump
- Flap designed to be opened and closed repeatedly
- Clear flap allows for better visibility
- Not made with natural rubber latex or plasticizer DEHP

Cat No.	Item	Qty/Unit
N705	Standard	40/box



Made in USA



# NeoBridge®

## Umbilical Catheter Holder



N700



N701



N702

### FEATURES AND BENEFITS

- Skin friendly NeoBond® Hydrocolloid base
- Secure positioning flaps designed to relieve stress on tubes and help prevent dislodgement
- Accommodates two catheters
- Flap designed to be opened and closed repeatedly
- Clear flap allows for better visibility
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Size	Approx Wt	Qty/Unit
N700	Micro	800 or less	20/box
N701	Small	800-1200 gm	20/box
N702	Large	1200-2200 gm	20/box



Made in USA U.S. Patent #5,370,627

# NeoSmile™

## Temperature Probe Cover



### FEATURES AND BENEFITS


- Skin friendly NeoBond® Hydrocolloid base
- Designed for long lasting adhesion
- Reflective foil helps protect against ambient and radiant heat
- Cute smiley face design
- Not made with natural rubber latex or plasticizer DEHP

Cat No.	Size	Qty/Unit
N731	1" Diameter	150/box



SoftSil™

Silicone Adhesive Patch



FEATURES AND BENEFITS

- Great for use as a skin barrier or guard
- Gentle silicone adhesive is ideal for neonatal skin
- Reposition as needed
- Not made with natural rubber latex or plasticizer DEHP

Cat No.	Item	Qty/Unit
N745	Silicone Adhesive Patch	100/box



# NeoPockets™


## Personal Incubator Organizer



### FEATURES AND BENEFITS

- Apply and remove repeatedly with no adhesive residue
- Crystal clear vinyl allows for visibility of patient
- Adheres to most smooth surfaces
- Three pockets for:
  - Patient info
  - Thermometer
  - Family photos and more
- Contours to a variety of incubator models
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Item	Qty/Unit
N7325	NeoPockets	20/box



# NeoSpoon™

## Disposable Measuring Spoon



Clinician Invented



1/8tsp (.6 ml)

1/2tsp (2.5 ml)



### FEATURES AND BENEFITS


- Dual spoon design (1/8 tsp and 1/2 tsp)
- Engineered for quick, accurate measuring
- Dispense exact volume in multiples of 1/8 tsp
- Ideal for feeding preparation
- Environmentally friendly, biodegradable material
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

Cat No.	Volume	Qty/Unit
N761	1/8 & 1/2 tsp	100/box



# Neotech Cell Wipe™

Cell Phone and Mobile Device Wipe



## FEATURES AND BENEFITS

- Disposable, pre-moistened wipe
- Lint free and non-abrasive
- Streak free and anti static
- Eliminates smudges, fingerprints and oils
- Ready to use dispenser box
- Individually packaged





# Neotech Cell Shield™

Cell Phone and Mobile Device Cover



## FEATURES AND BENEFITS

- Helps to keep cell phones and other mobile devices clean
- Touch, talk, swipe and take pictures right through the clear cover
- Disposable
- Ready to use dispenser box
- Not made with natural rubber latex or plasticizer DEHP

Cat No.	Item	Qty/Unit
N15250	Cell Wipe	1000/case*
**4 boxes of 250		
N15100	ComboPack with Large Shield	400ea/case**
***4 boxes of 100 Cell Shields & 100 Cell Wipes		



# NeoGrip®

## Tubing and Cable Holder



### FEATURES AND BENEFITS

- Helps to keep tubes, cables and circuits together
- Allows for better patient mobility
- Attaches to linen or clothing without damaging
- Helpful during kangaroo care
- Open and close repeatedly with no residue
- Can be used in clinical or home care setting
- Great for ambulatory patients
- Not made with natural rubber latex or plasticizer DEHP

U.S. Patent #D814,568S



Cat No.	Color	Qty/Unit
N601	Teal (VersaGrip)	40/box
N602	Pink/Blue	40/box
N602P	Pink	40/box
N602B	Blue	40/box
N602R	Red	40/box
N602G	Green	40/box
N602Y	Yellow	40/box
N602O	Orange	40/box
N602V	Violet	40/box
N602A	Assorted	40/box
N603	White	40/Box

## Clinical Training

Neotech employs an experienced team of licensed Clinical Consultants who work with hospitals providing educational support and in-service training.

Our highly specialized team is located throughout the United States. Their expertise in the neonatal, pediatric and respiratory fields is exercised to educate you on the safe and proper use of our products.




Clinical Consultants can provide onsite in-service training for your staff including night and weekend shifts. Our Clinical Consultants are also available to assist our distributor reps.

To schedule your on-site training, call **1-800-966-0500** or visit **[neotechproducts.com/clinical-training](https://neotechproducts.com/clinical-training)**.

## Neotech Research

Neotech Products offers two unique opportunities to support clinicians, like you, who are conducting neonatal and pediatric clinical research.

The **Professional Development Grant Program** awards monetary assistance to clinicians so that they can present their project results as podium or poster sessions at local, national, and international meetings.



Apply for the **Neotech Research Support Grant** to receive monetary and/or product support for unit-based clinical evaluations, quality improvement projects using any Neotech product; as well as IRB-approved research projects of any size.

Visit [neotechresearch.com](https://neotechresearch.com) to learn more and apply.



**For more information or to order, call**

*our Sole Distributor in Ireland*

**Pat Coakley Medical Ltd**

Units E3 E4  
Network Enterprise Park  
Kilcoole A63 H244  
Co Wicklow  
IRELAND

**+353 1 2018070**  
**sales@patcoakleymedical.ie**

